

# EXTENSION PACK FOR GAMMA PLATFORM

2025.1

|                                   |    |
|-----------------------------------|----|
| <b>Hardware requirements</b>      | 3  |
| What is needed?                   | 3  |
| <b>Therapeutic tasks database</b> | 5  |
| Speed                             | 5  |
| Balance                           | 9  |
| Movement precision                | 11 |
| Functional movements              | 13 |

# WHAT IS NEEDED?

Please make sure the PC where you want this module to be active have VAST.Rehab Patient Panel installed and that the following hardware requirements are met:

- Windows 10/11
- Intel Core i5 (8th gen or newer). Important: Avoid ultra-low-power versions (e.g., i5-8250U), as they may not meet performance requirements. Prefer standard or high-performance CPUs.
- Minimum: 8 GB RAM (16 GB or more recommended for optimal performance).




# SPEED


## RABBIT

Measure number of repetitions of specific movement pattern an individual is able to perform within predefined time interval.

### CONTROL MODES



### RESULTS



### ADJUSTMENTS

- Task duration
- Height of jump

### OBJECTIVES

- Speed of movement
- Repetitive movements

### INSTRUCTION FOR PATIENT

Go through the entire route as fast as you can.





## SPEED

RABBIT

### SAMPLE SETTINGS





Duration

< 90s >

Height of jump

< 10 >




# SPEED


## KAYAK

Measure number of repetitions of specific movement pattern an individual is able to perform within predefined time interval.

### CONTROL MODES



### RESULTS



### ADJUSTMENTS

- Task duration
- Height of jump

### OBJECTIVES

- Speed of movement
- Repetitive movements

### INSTRUCTION FOR PATIENT


Row as fast as you can.



## SAMPLE SETTINGS



|          |                |
|----------|----------------|
|          |                |
| Duration | Height of jump |
| < 90s >  | < 10 >         |




# BALANCE

## BLOCK BUILDER

Measure and train individual's skills to perform specific movement patterns while keeping predefined weight distribution.

## CONTROL MODES



## RESULTS



## ADJUSTMENTS

- Task duration
- Height of jump
- Stack height

## OBJECTIVES

- Movement precision
- Muscle strengthening


## INSTRUCTION FOR PATIENT

Build as many stacks as you can. Keep your body balanced.





## SAMPLE SETTINGS



|                 |                          |                      |
|-----------------|--------------------------|----------------------|
|                 |                          |                      |
|                 | Difficulty<br><b>1/3</b> |                      |
| Duration<br>90s |                          | Height of jump<br>10 |
|                 |                          | Stack height<br>6    |



# MOVEMENT PRECISION


## UMBRELLA

Measure and train individual's skills to perform specific movement patterns with predefined speed and range.

## CONTROL MODES



## RESULTS



Total points

14 point(s)

Movement precision

32 %

## ADJUSTMENTS

- Task duration
- Path
- Range
- Umbrella size

## OBJECTIVES

- Movement precision
- Visual motor coordination

## INSTRUCTION FOR PATIENT

Don't let the hippo get wet - keep the umbrella above him!



# MOVEMENT PRECISION

## UMBRELLA

### SAMPLE SETTINGS



|                        |                          |              |
|------------------------|--------------------------|--------------|
|                        |                          |              |
|                        | Difficulty<br><b>1/3</b> |              |
| Duration<br><b>60s</b> |                          | Path<br>8.0s |
| Range<br>0% ↔ 100%     | Umbrella size<br>150%    |              |




# FUNCTIONAL MOVEMENTS


## DRAGON

Measure and train individual's skills to perform movements based on real-world situational biomechanics. They usually involve multi-planar, multi-joint movements which place demand on the body's core musculature and innervation.

## CONTROL MODES



## RESULTS



Total points



20 point(s)



Functional movements



23 %

## ADJUSTMENTS

- Task duration
- Coins group size
- Height of jump
- Distance between coins
- Gravity force

## OBJECTIVES

- Predicting the trajectory of objects
- Improve range of motion
- Visual motor coordination
- Muscle strengthening
- Planning and Strategy

## INSTRUCTION FOR PATIENT

Fly and collect the coins.




# FUNCTIONAL MOVEMENTS

## DRAGON

### SAMPLE SETTINGS



|                                    |                             |                           |
|------------------------------------|-----------------------------|---------------------------|
|                                    |                             |                           |
| ◀                                  | Difficulty<br><b>custom</b> | ▶                         |
| Duration<br>< 90s >                |                             | Coins group size<br>< 3 > |
| Distance between coins<br>< 250% > |                             | Height of jump<br>< 10 >  |
| Gravity force<br>< 100% >          |                             |                           |



|                                    |                          |                           |
|------------------------------------|--------------------------|---------------------------|
|                                    |                          |                           |
| ◀                                  | Difficulty<br><b>1/3</b> | ▶                         |
| Duration<br>< 90s >                |                          | Coins group size<br>< 5 > |
| Distance between coins<br>< 250% > |                          | Height of jump<br>< 10 >  |
| Gravity force<br>< 100% >          |                          |                           |




# FUNCTIONAL MOVEMENTS


## RAILS

Measure and train individual's skills to perform movements based on real-world situational biomechanics. They usually involve multi-planar, multi-joint movements which place demand on the body's core musculature and innervation.

### CONTROL MODES



### RESULTS



### ADJUSTMENTS

- Speed
- Task duration
- Range
- Route shape
- Enable derailing
- Enable obstacles
- Time between objects

### OBJECTIVES

- Dynamic responses to emerging moving targets
- Predicting the trajectory of objects
- Visual motor coordination

### INSTRUCTION FOR PATIENT

Control the trolley to collect the coins.





# FUNCTIONAL MOVEMENTS

## RAILS

### SAMPLE SETTINGS





◀

Difficulty

▶

1/3

Speed

< 100% >

speed set automatically

Duration

< 90s >

Range

< 0% ↔ 100% >

Route shape

< >

Enable derailling


< No >



Enable obstacles

< No >

Time between objects

< 5s >





◀

Difficulty

▶

custom

Speed

< 100% >

speed set automatically

Duration

< 90s >

Range

< 0% ↔ 100% >

Route shape

< >

Enable derailling

< No >

Enable obstacles

< No >

Time between objects

< 5s >



# FUNCTIONAL MOVEMENTS


## MOTOCROSS

Measure and train individual's skills to perform movements based on real-world situational biomechanics. They usually involve multi-planar, multi-joint movements which place demand on the body's core musculature and innervation.

## CONTROL MODES



## RESULTS



Total points

29 point(s)



Functional movements

14 %

## ADJUSTMENTS

- Task duration
- Range
- Route shape

## OBJECTIVES

- Dynamics of planned movements
- Planning and Strategy

## INSTRUCTION FOR PATIENT

Accelerate and brake to cover the entire route as quickly as possible without tipping.





## SAMPLE SETTINGS



|                     |                          |                    |
|---------------------|--------------------------|--------------------|
|                     |                          |                    |
|                     | Difficulty<br><b>1/3</b> |                    |
| Duration<br>90s     |                          | Range<br>0% ↔ 100% |
| Route shape<br>Easy |                          |                    |



# FUNCTIONAL MOVEMENTS


## FOREST RUNNER

Measure and train individual's skills to perform movements based on real-world situational biomechanics. They usually involve multi-planar, multi-joint movements which place demand on the body's core musculature and innervation.

## CONTROL MODES



## RESULTS



## ADJUSTMENTS

- Speed
- Task duration
- Range

## OBJECTIVES

- Dynamics of planned movements
- Focusing
- Planned movements
- Speed of movement

## INSTRUCTION FOR PATIENT

Keep the hare on the run, avoid obstacles and collect as many carrots as you can.





# FUNCTIONAL MOVEMENTS

## FOREST RUNNER

### SAMPLE SETTINGS





◀

Difficulty  
**1/2**

▶

<

Speed  
**150%**

>


speed set automatically

<


Duration  
**90s**

>

<

Range  
  
0% ↔ 100%

>



# FUNCTIONAL MOVEMENTS


## DANCEMAN

Measure and train individual's skills to perform movements based on real-world situational biomechanics. They usually involve multi-planar, multi-joint movements which place demand on the body's core musculature and innervation.

## CONTROL MODES



## RESULTS



Total points

27 point(s)



Functional movements

12 %

## ADJUSTMENTS

- Task duration
- Range
- Advanced scoring
- Song index
- Spawn rate level

## OBJECTIVES

- Activity in a given rhythm
- Spontaneous movements
- Visual motor coordination

## INSTRUCTION FOR PATIENT

Hit the green characters when they come close.



# FUNCTIONAL MOVEMENTS

## DANCEMAN

### SAMPLE SETTINGS



|                                 |                          |
|---------------------------------|--------------------------|
|                                 |                          |
|                                 | Difficulty<br><b>1/6</b> |
| Duration<br><b>90s</b>          | Range<br><br>0% ↔ 100%   |
| Advanced scoring<br><b>No</b>   | Song index<br><b>0</b>   |
| Spawn rate level<br><b>Easy</b> |                          |