

RELAXATION PACK

2025.1

Hardware requirements	3
What is needed?	3
Therapeutic tasks database	5
Functional movements	5

WHAT IS NEEDED?

Please make sure the PC where you want this module to be active have VAST.Rehab Patient Panel installed and that the following hardware requirements are met:

- Oculus Quest 2
- Oculus Quest 3




FUNCTIONAL MOVEMENTS


SEA CREATURES

Measure and train individual's skills to perform movements based on real-world situational biomechanics. They usually involve multi-planar, multi-joint movements which place demand on the body's core musculature and innervation.

CONTROL MODES



RESULTS



ADJUSTMENTS

- Task duration

OBJECTIVES

INSTRUCTION FOR PATIENT

Feel the experience of swimming in the deep sea surrounded by friendly sea creatures. Use primary index trigger button (or A/X button) to capture a photo to make creatures react and swim closer to you. Aim at the center of the creature for a higher photo score.




FUNCTIONAL MOVEMENTS


ARCHEOLOGY

Measure and train individual's skills to perform movements based on real-world situational biomechanics. They usually involve multi-planar, multi-joint movements which place demand on the body's core musculature and innervation.

CONTROL MODES



RESULTS



ADJUSTMENTS

- Task duration
- Stone hardness
- Hand

OBJECTIVES

INSTRUCTION FOR PATIENT

Position the chisel atop the stone, then strike it using the hammer. This action will break apart the stone, unveiling the hidden artifact within.




FUNCTIONAL MOVEMENTS


PLANTING TREES

Measure and train individual's skills to perform movements based on real-world situational biomechanics. They usually involve multi-planar, multi-joint movements which place demand on the body's core musculature and innervation.

CONTROL MODES



RESULTS



ADJUSTMENTS

- Task duration
- Hand

OBJECTIVES

INSTRUCTION FOR PATIENT

Launch seeds using a slingshot at the targets highlighted by flying objects. Use your controller's joystick to choose from a variety of tree colors.